


1. Medical mask gas exchange pressure difference tester



Use:

This tester is used to determine the ventilation resistance and pressure difference indicators of medical surgical masks under specified conditions.

Features:

1. HD color touch screen.
2. High precision imported digital differential pressure meter.
3. High precision imported flowmeter, high precision flow control.

Parameters:

Flowmeter range	0-10L/min, Precision:3%
Digital differential pressure measurement range	0-2000Pa, Precision:1Pa
Air volume of the air pump	8L/min
Experimental area	4.9cm ² 、5.06cm ²
External dimension	460×500×670mm(L×W×H)
Weight	25kg
power supply	AC 220V 50Hz 300W

Standards:

YY 0469-2004

Y/T 0969-2013